Cleaning Instructions

Foam Overlay

Instructions de nettoyage du

Superposition de mousse

Blake
Medical Distribution

1-866-662-3459
www.blakemedical.ca

For daily cleanings all that is needed is a wipe down of the cover. This cover is a fluid proof fabric and will prevent most fluids from permeating through to the overlay. For these daily cleanings soap and water will be all that is needed. If fluids do soil the overlay follow the instructions below.

Covers

Remove all contamination with cloth or paper towel, then wipe the surface with any of the following. 3M Hard Surface Disinfectant Cleaner Wipes, Wood Wyant ED Disinfectant, Virox 5, Oxvir TB, Virex ® II 256 One-step Disinfectant Cleaner & Deodorant, Virex II 128 One-step Disinfectant Cleaner, Sanikleen and Super Sani Cloth Germicidal Disposable Wipe. Next disinfect liberally with any of the above wipes (avoid pooling of the disinfectant), wipe dry after 1-3 minutes. Next wipe the surface with a damp cloth and mild soap and finally a damp cloth with liberal amounts of water. The cover is machine washable, up to 60 degrees C. Do not wring, lay flat to dry.

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Gel Infused Foam Overlay Cleaning

Our Gel Infused Foam overlay is not washable ensure that the foam does not get wet. The incontinent cover will protect the foam from becoming wet, but if the foam overlay becomes heavily soiled, a replacement will need to be purchased.

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Care of your Overlay

Blake Medical's Overlays are designed to be used at room temperature. If it is subjected to extreme heat or cold, allow the overlay to acclimate to room temperature before using.

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