


Products Cushion Comparison Chart

	GEL					FOAM		FOAM & GEL							
	BRONZE	SILVER	GOLD	GOLD+	BOLD	ZENTURY	PC Cutback	HYBRID	HYBRID LP	HYBRID DEEP	HY. DEEP LP	SPLIT	PLATINUM	PLAT. PLUS	
Great For:															
Risk of Skin Breakdown	★	★★★★	★★★★	★★★★	★★★★	★	★	★★★★	★★	★★★★	★★	★★★★	★★★★	★★★★	
Shape/Contour	★	★	★	★	★	★★★★	★★	★★	★★	★★★★	★★★★	★★	★★	★★	
Immersion	★★	★★★★	★★★★	★★★★	★★★★	★★★★	★★★★	★★	★★	★★★★	★★	★★	★★★★	★★★★	
Envelopment	★★	★★★★	★★★★	★★★★	★★★★	★★★★	★★	★★	★★	★★★★	★★★★	★★	★★★★	★★★★	
Stability	★★	★★	★★	★★	★★	★★	★★	★★★★	★★★★	★★★★	★★★★	★★★★	★★	★★	
Shear reduction	★★★★	★★★★	★★★★	★★★★	★★★★	★	★	★★	★★	★★	★★	★★	★★★★	★★★★	
Vibration Dampening	★★	★★★★	★★★★	★★★★	★★★★	★	★	★★	★★	★★	★★	★★	★★★★	★★★★	
Adjustment	★★★★	★★★★	★★★★	★★★★	★★★★	★		★★	★★	★★	★★	★★★★	★★	★★	
Temp Control	★★	★★	★★	★★	★★	★★★★	★★★★	★	★	★	★	★	★	★	
Product Weight	★★	★★★★	★★★★	★★★★	★★★★	★	★	★★	★★	★★	★★	★★	★★★★	★★★★	
Weight Capacity	★	★	★	★★	★★★★	★	★	★	★	★	★	★	★★	★★★★	


Products

Cushion Comparison Chart


High Risk




200lbs
Geo-Matrix Silver




250lbs
Geo-Matrix Gold



300lbs
Geo-Matrix Gold+




600lbs
Geo-Matrix Bold




450lbs
Geo-Matrix Platinum+


Multi-user – Suitable for Moderate & High Risk




275lbs
Hybrid



275lbs
Hybrid Deep




275lbs
Geo-Matrix Platinum




125kg
Geo-Matrix Split

Moderate Risk




250lbs Geo-Matrix Hybrid
Deep Low Profile




250lbs Geo-Matrix
Hybrid Low Profile


Low Risk



275lbs
Zentury



275lbs
PC Cutback



200lbs
Geo-Matrix Bronze

Blake Medical realizes that everyone is different and requires different cushion attributes to make each individual comfortable and safe. Please note that the above chart is for reference only and that each individual user should contact a healthcare professional to ensure which cushion is right for them.